

Orthopaedic Studio v1.2.1 User Manual

<http://orthostudio.spectronic.se>

1. Orthopaedic Studio - Instructions for usage	3
1.1 Introduction	3
1.2 Installation	3
1.3 General Usage	3
1.4 Update log	4
2. Description of measured parameters	5
2.1 Anteroposterior visual scores.....	5
2.1.1 Break in Shenton's line	5
2.1.2 Cross-over sign.....	5
2.1.3 Posterior wall sign.....	6
2.1.4 Tönnis classification.....	6
2.2 Anteroposterior quantitative measurements.....	6
2.2.1 LCE (Lateral Center Edge) angle	6
2.2.2 Tönnis angle	7
2.2.3 Joint space width	7
2.2.4 Pelvic tilt and rotation.....	7
2.3 Von Rosen visual scores	8
2.3.1 Joint Congruity.....	8
2.4 False profile quantitative measurements.....	8
2.4.1 ACE (Anterior Center Edge) angle.....	8
2.5 Frog image quantitative measurements.....	9
2.5.1 Alpha angle	9
2.5.2 Alpha angle for SCFE patients	9
2.5.3 Head-neck offset ratio.....	9
2.5.4 Southwick angle.....	10
2.5.5 Epiphysiseal-metaphyseal offset	10
3. Disclaimer	11

1. Orthopaedic Studio - Instructions for usage

1.1 Introduction

Orthopaedic Studio is a plugin for OsiriX, designed to help orthopaedic specialists perform a number of common quantitative hip examinations based on standard radiograph images.

Orthopaedic Studio is designed for ease of use. The user identifies landmarks in the radiograph images, according to the self-explanatory instructions and illustrations which are shown at each step of the process. From these landmarks Orthopaedic Studio automatically calculates all angles, measurement and offsets that are of interest in the image. The final results can be directly exported into Excel compatible files.

1.2 Installation

Orthopaedic Studio can be downloaded and installed using the built-in 'Plugins manager' in OsiriX, which is found in the Plugins menu. Click on the downloads tab and then select to download and install Orthopaedic Studio.

Orthopaedic Studio can also be downloaded from our website (<http://orthostudio.spectronic.se/>).

1.3 General Usage

Orthopaedic Studio can do hip evaluations on four different types of radiograph images: Anteroposterior (AP), Von Rosen, False Profile and Frog.

The procedure is self-explanatory and very similar for all types of images. The example below describes how to evaluate an AP image:

1. First open the AP image to be evaluated in OsiriX.
2. Start Orthopaedic Studio by selecting 'Anteroposterior' from the Orthopaedic Studio sub menu. (*Plugins -> ROI Tools -> Orthopaedic Studio -> Anteroposterior*).
3. Select the appropriate visual scores for the current image. Then click the 'Set markers' button.
4. Place markers to identify landmarks within the image, in accordance with the on-screen instructions. Note that the markers can be moved after being placed. When all markers are placed, click the 'Results' button.
5. Review the results and perform changes if necessary.
6. Click on any of the save buttons to save the data to an Excel compatible file.

The other types of evaluations (i.e. Von Rosen, False Profile or Frog) are performed according to similar schemes. Substitute 'AP' for any of the other evaluation types in steps 1 and 2 above and follow the on-screen instructions.

Note that Orthopaedic Studio exports its Excel data into a '.csv' file instead of an '.xls' file. The '.csv' file is fully compatible with Excel and can be opened by double-clicking the file in a Finder window.

1.4 Update log

Version 1.0.0

- Initial release

Version 1.1.0

- Added short keys for starting the plugin (Shift-Ctrl-A for AP, Shift-Ctrl-F for Frog, Shift-Ctrl-V for Von Rosen, Shift-Ctrl-P for False Profile).
- No longer necessary to manually close the plugin window before changing between evaluation modes.
- Save file has been reorganized. Every evaluation is now placed on a new line.
- Labels for markers are not shown while they are initially being placed.
- ACE angle can now be negative.
- Tönnis angle can now be negative.
- Instructions for placing marker on pubic symphysis has been updated.

Version 1.2.0

- Added support for evaluation of SCFE patients
- License manager included in plugin
- Several minor improvements

Version 1.2.1

- Added support for Mac OS X 10.7
- Added a help button from within the plugin


2. Description of measured parameters

2.1 Anteroposterior visual scores

2.1.1 Break in Shenton's line

There is an imaginary line connecting the medial aspect of the femoral neck to the superior pubic ramus. It is used to determine the relationship/alignment of the head of the femur to the acetabulum.


Observe the alignment of the femur and pelvis and determine if they align in an approximately smooth arc. If the curve is "broken" by more than 5mm then Shenton's line is said to have a break. Interruption of the curve occurs with subluxation or dislocation of the hip.



Reference: Siebenrock, K.A., D.F. Kalbermatten, and R. Ganz, *Effect of pelvic tilt on acetabular retroversion: a study of pelves from cadavers*. Clin Orthop Relat Res, 2003(407): p. 241-8.

2.1.2 Cross-over sign

A hip with a normal pelvic inclination should have the anterior and posterior lips join at the edge of the acetabulum. If a cross is formed by the radiographic projections of the anterior and posterior lips then the hip is positive for "cross-over sign".




Reference: Reynolds, D., J. Lucas, and K. Klaue, *Retroversion of the acetabulum. A cause of hip pain*. J Bone Joint Surg Br, 1999. 81(2): p. 281-8.

2.1.3 Posterior wall sign

If the posterior wall of acetabulum is medial to the center of the femoral head then hip is positive to "posterior wall sign"

Reference: Reynolds, D., J. Lucas, and K. Klaue, *Retroversion of the acetabulum. A cause of hip pain.* J Bone Joint Surg Br, 1999. 81(2): p. 281-8.



2.1.4 Tönnis classification

A subjective categorization of the presence of osteoarthritis in the hip joint. The classification involves four categories with the following criterion:

- Grade 0: Normal, no degenerative changes, no signs of OA
- Grade 1: Subchondral Sclerosis, minimal JS narrowing, and minimum formation of osteophytes
- Grade 2: Hips with subchondral cyst formation and moderate JS narrowing
- Grade 3: Hips with severe or complete but localized JS narrowing
- Grade 4: Hips that have undergone extensive or severe cartilage loss.


Reference: Clohisy, J.C., et al., *A systematic approach to the plain radiographic evaluation of the young adult hip.* J Bone Joint Surg Am, 2008. 90 Suppl 4: p. 47-66.

2.2 Anteroposterior quantitative measurements

2.2.1 LCE (Lateral Center Edge) angle

The angle formed between the two lines passing through the center of the femoral head, one of which extends to the lateral edge of sourcil (A) and the line perpendicular (B) to that joining the centers of the two femoral heads (of the two hips). The normal angle of an adult is greater than 25°. Beyond 40° may indicate pincer impingement.


Reference: Wiberg, G., *Studies on dysplastic acetabula and congenital subluxation of the hip joint.* Acta Chir Scand, 1939(58 (suppl)): p. 5-135.



2.2.2 Tönnis angle

The angle formed by the intersection of horizontal line connecting the femoral head centers and the line that passes through medial edge of the sourcil. Normal is 0° - 10° , less than 0° indicates risk of impingement and higher than 10° indicates structural instability.

Reference: Clohisy, J.C., et al., *A systematic approach to the plain radiographic evaluation of the young adult hip*. J Bone Joint Surg Am, 2008. 90 Suppl 4: p. 47-66.



2.2.3 Joint space width

The smallest distance between the femoral head and the acetabulum. The distance should normally be 3-5mm. Smaller distance than 3mm, or a distance that is smaller than that of the opposite hip, indicates radiological osteoarthritis.

Reference: Fredensborg N, Nilsson BE. *The joint space in normal hip radiographs*. Radiology 1978;126:325-6

2.2.4 Pelvic tilt and rotation

Distance between pubic symphysis and the tip of coccyx (i.e. pelvic tilt) should be 1–3 cm. Tip of coccyx should be in line with pubic symphysis (i.e. pelvic rotation should be low).

Reference: Siebenrock, K.A., D.F. Kalbermatten, and R. Ganz, *Effect of pelvic tilt on acetabular retroversion: a study of pelves from cadavers*. Clin Orthop Relat Res, 2003(407): p. 241-8.

2.3 Von Rosen visual scores

2.3.1 Joint Congruity

Visual classification of joint congruency. Can be used both with images taken in the Von Rosen position and standard standing AP position.


Excellent - if the curvature of the acetabulum and that of the femoral head are almost identical and the joint space is adequately maintained.

Good - if the curvature of the acetabulum and the femoral head are not identical but the joint space is adequately maintained.

Fair - if partial narrowing of the joint space has occurred.

Poor - if partial disappearance of joint space has occurred.

Reference: Yasunaga, Y., et al., *Rotational acetabular osteotomy in patients forty-six years of age or older: comparison with younger patients*. J Bone Joint Surg Am, 2003. 85-A(2): p. 266-72.




2.4 False profile quantitative measurements

2.4.1 ACE (Anterior Center Edge) angle

A measure to determine the anterior coverage of the femoral head and to diagnose arthritis affecting the anterior part of the joint.

The ACE angle is measured on images in the false profile projection and is defined as the angle between an imagined vertical line and the line that extends from the center of the femoral head to the lateral edge of the sourcil. An ACE angle less than 20° is indicative for structural instability.

Reference: Lequesne, M.G. and J.D. Laredo, *The faux profil (oblique view) of the hip in the standing position. Contribution to the evaluation of osteoarthritis of the adult hip*. Ann Rheum Dis, 1998. 57(11): p. 676-81.




2.5 Frog image quantitative measurements

2.5.1 Alpha angle

The alpha angle is the angle between two lines originating in the center of the femoral head. One line extends out through the center of the femoral neck and the other line extends out through the point where the radius of the head starts to increase compared the central aspect (i.e. where the femoral head extends out of its boundary circle).

An alpha angle beyond 55° indicates cam impingement.


Reference: Notzli, H.P., et al., *The contour of the femoral head-neck junction as a predictor for the risk of anterior impingement*. J Bone Joint Surg Br, 2002. 84(4): p. 556-60.



2.5.2 Alpha angle for SCFE patients

Due to the translation of the femoral head in SCFE (slipped capital femoral epiphysis) patients, a modified alpha angle is sometimes preferred instead of the ordinary alpha angle. The modification means that the line drawn through the femoral neck does not have to also go through the center of the femoral head. Instead, this line should just be set parallel to the femoral neck.


Reference: Leunig, M., et al, *In Situ Pinning With Arthroscopic Osteoplasty for Mild SCFE*, Clin Orthop Relat Res, 2010. 468: p.3160–67



2.5.3 Head-neck offset ratio

The head-neck offset ratio is calculated as the distance between lines A and B, divided by the diameter of the femoral head, C. Line A extends through the anteriormost aspect of the femoral neck, while line B extends through the anteriormost aspect of the femoral head. Both lines are parallel to the long axis of the femoral neck. An offset ratio that is less than 0.17 indicates that cam deformity is likely.


Reference: Peelle, M.W., et al., *Acetabular and femoral radiographic abnormalities associated with labral tears*. Clin Orthop Relat Res, 2005. 441: p. 327-33.



2.5.4 Southwick angle

The Southwick angle is of interest in determining the amount of correction required for SCFE (slipped capital femoral epiphysis) cases. It is defined as the angle between lines A and B, where line A extends in the direction of the femoral shaft and line B is perpendicular to line C, which connects the edges of the slipped epiphysis.


Reference: Southwick, W.O., *Osteotomy through the Lesser Trochanter for Slipped Capital Femoral Epiphysis*, J Bone Joint Surg Am, 1967. 49: p.807-35.



2.5.5 Epiphysiseal-metaphyseal offset

The epiphysiseal-metaphyseal offset is measured as the distance between a line drawn along the anterior edge of the epiphysis parallel to the femoral neck and a line drawn along the proximal-anterior edge of the metaphysis also parallel to the femoral neck.

Reference: Leunig, M., et al, *In Situ Pinning With Arthroscopic Osteoplasty for Mild SCFE*, Clin Orthop Relat Res, 2010. 468: p.3160–67



3. Disclaimer

The Orthopaedic Studio software and its documentation is provided as is. Use the software and its documentation at your own risk. The software does not have an FDA clearance and is therefore restricted to investigational use.

NEITHER THE AUTHOR NOR SPECTRONIC MAKE ANY WARRANTIES AS TO PERFORMANCE OR FITNESS FOR A PARTICULAR PURPOSE, OR ANY OTHER WARRANTIES WHETHER EXPRESSED OR IMPLIED. NO ORAL OR WRITTEN COMMUNICATION FROM OR INFORMATION PROVIDED BY THE AUTHOR OR SPECTRONIC SHALL CREATE A WARRANTY. UNDER NO CIRCUMSTANCES SHALL THE AUTHOR OR SPECTRONIC BE LIABLE FOR DIRECT, INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES RESULTING FROM THE USE, MISUSE, OR INABILITY TO USE THE SOFTWARE, EVEN IF THE AUTHOR OR SPECTRONIC HAVE BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

ORTHOPAEDIC STUDIO IS NOT INTENDED FOR PRIMARY DIAGNOSTICS.

All calculations, measurements and images provided by this software are intended for investigational use only.

If you use Orthopaedic Studio for scientific research please give appropriate credit in publications.

<http://orthostudio.spectronic.se>